

# Syntaxin 10 antibody



Catalog Number: 115790

## Product name

Syntaxin 10 antibody

## Immunogen

[Human STX10 Recombinant protein \(GST tag & His tag\)](#)

## Specificity

Human; other species not tested.

## Antibody description

Syntaxin 10 Rabbit Polyclonal antibody. Positive WB detected in BxPC-3 cells. Observed molecular weight by Western-blot: 28 kDa

## Preparation

This antibody was obtained by immunization of Syntaxin 10 recombinant protein (Accession Number: NM\_001271611). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:200-1:1000

## Validations



BxPC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:115790(STX10 antibody) at dilution of 1:200