

Anti-DDX3 (6G8) Mouse antibody



Catalog Number: 168408

Product name

Anti-DDX3 (6G8) Mouse antibody

Specificity

Human, Rat, Mouse, Monkey

Antibody description

DDX3 (6G8) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant human DDX3 protein fragments expressed in E.coli.

Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.02% sodium azide 0.1mg/mlBSA and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG2a

Applications


WB, IP, ICC/IF, ChIP

Dilutions


WB: 1/1000

;ICC: 1/200

Validations




Western blot detection of DDX3 in HeLa, 3T3, C6, COS7, K562 and Jurkat cell lysate using DDX3 mouse mAb (1:1000 diluted). Predicted band size: 75KDa. Observed



band size: 75KDa.

Immunocytochemistry staining of HeLa cells fixed with 4% Paraformaldehyde and using DDX3



mouse mAb (dilution 1:200).

Immunoprecipitation analysis of HeLa cell lysates using DDX3 mouse mAb.