

TMEM176B antibody



Catalog Number: 116108

Product name

TMEM176B antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

TMEM176B Rabbit Polyclonal antibody. Positive IP detected in mouse lung tissue. Positive WB detected in A549 cells, human placenta tissue. Observed molecular weight by Western-blot: 31 kDa, 23 kDa

Preparation

This antibody was obtained by immunization of TMEM176B recombinant protein (Accession Number: XM_047420259). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IP: 1:200-1:2000

Validations



A549 cells were subjected to SDS PAGE followed by western blot with Catalog No:116108(TMEM176B Antibody) at dilution of



1:300

IP Result of anti-TMEM176B (IP:Catalog No:116108, 4ug; Detection:Catalog No:116108 1:500) with mouse lung tissue lysate 4000ug.