

# RAB8A Polyclonal Antibody



Catalog Number: 163508

---

## Product name

RAB8A Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to RAB8A

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-207 of human RAB8A (NP\_005361.2).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

IF

## Dilutions

IF 1:50 - 1:200

## Validations



Immunofluorescence - RAB8A Polyclonal Antibody

Immunofluorescence analysis of U-2 OS cells using RAB8A antibody at dilution of 1:100. Blue: DAPI for nuclear staining.