

TESK2 Polyclonal Antibody



Catalog Number: 161608

Product name

TESK2 Polyclonal Antibody

WB 1:500 - 1:2000

IF 1:50 - 1:200

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to TSK2

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 302-571 of human TSK2 (NP_009101.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type


Rabbit IgG

Applications

WB, IF


Dilutions

Validations



Western blot - TSK2 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using TSK2 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 1s.



Immunofluorescence - TSK2 Polyclonal Antibody

Immunofluorescence analysis of A549 cells using TSK2 antibody. Blue: DAPI for nuclear staining.