

TRIM27 antibody



Catalog Number: 116308

Product name

TRIM27 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TRIM27 Rabbit Polyclonal antibody. Positive IHC detected in human lung cancer tissue, human brain tissue, human heart tissue, human kidney tissue, human lung tissue, human ovary tissue, human placenta tissue, human skin tissue, human spleen tissue, human testis tissue. Positive IF detected in HeLa cells. Positive WB detected in Jurkat cells, human testis tissue, mouse testis tissue. Positive IP detected in mouse testis tissue. Observed molecular weight by Western-blot: 58 kDa

Preparation

This antibody was obtained by immunization of TRIM27 recombinant protein (Accession Number: NM_006510). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC, IP, IF

Dilutions

Recommended Dilution:


WB: 1:500-1:2000

IP: 1:500-1:5000


IHC: 1:20-1:200

IF: 1:20-1:200

Validations




Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:116308(TRIM27 antibody) at dilution of 1:1000




Immunohistochemical of paraffin-embedded human lung cancer using Catalog No:116308(TRIM27 antibody) at dilution of 1:100 (under 10x lens)

TRIM27 antibody

Catalog Number: 116308



IP Result of anti-TRIM27 (IP:Catalog No:116308, 4ug; Detection:Catalog No:116308 1:1000) with mouse testis tissue lysate 4000ug.



Immunofluorescent analysis of HeLa cells using Catalog No:116308(TRIM27 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)