

FOXP1 antibody



Catalog Number: 110708

Product name

FOXP1 antibody

Immunogen

[Human FOXP1 Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

FOXP1 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, A549 cells, LNCaP cells, MCF7 cells, MDA-MB-453s cells, mouse heart tissue, mouse lung tissue, mouse testis tissue, PC-3 cells, Raji cells, Y79 cells. Observed molecular weight by Western-blot: 50kDa, 60-65 kDa, 85kDa

Preparation

This antibody was obtained by immunization of FOXP1 recombinant protein (Accession Number: BC131720). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:110708(FOXP1 antibody) at dilution of 1:1000