

VAMP5 antibody



Catalog Number: 116708

Product name

VAMP5 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

VAMP5 Rabbit Polyclonal antibody. Positive WB detected in human lung tissue, mouse kidney tissue, PBMC cells. Positive IHC detected in human lung cancer tissue. Observed molecular weight by Western-blot: 13-15kd

Preparation

This antibody was obtained by immunization of VAMP5 recombinant protein (Accession Number: NM_006634). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



human lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:116708(VAMP5 antibody) at dilution of 1:400



WB result of VAMP5 antibody (Catalog No:116708, 1:300) with PBMC cells. Proteintech α VAMP5 in 10% BSA / TBS 0.05% Tween-20 O/N @ 4°C Goat anti-rabbit HRP 1:10000 in 5% milk / TBS 0.05% Tween-20 1h @ RT.



Immunohistochemical of paraffin-embedded

VAMP5 antibody



Catalog Number: 116708

human lung cancer using Catalog

No:116708(VAMP5 antibody) at dilution of 1:50
(under 10x lens)