

SCRG1 antibody



Catalog Number: 115018

Product name

SCRG1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

SCRG1 Rabbit Polyclonal antibody. Positive IHC detected in mouse brain tissue, human testis tissue.

Preparation

This antibody was obtained by immunization of SCRG1 recombinant protein (Accession Number: BC017583). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200


Validations



Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:115018(SCRG1 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:115018(SCRG1 Antibody) at dilution of 1:50 (under 40x lens)




Immunohistochemistry of paraffin-embedded mouse brain tissue slide using Catalog No:115018(SCRG1 Antibody) at dilution of 1:50 (under 10x lens)

SCRG1 antibody



Catalog Number: 115018



Immunohistochemistry of paraffin-embedded mouse brain tissue slide using Catalog No:115018(SCRG1 Antibody) at dilution of 1:50 (under 40x lens)