

Recombinant CHMP5 protein (His tag)



Catalog Number: 541418

Gene Name Synonym

Chromatin-modifying protein 5; SNF7 domain-containing protein 2; Vacuolar protein sorting-associated protein 60; Vps60; hVps60

Protein Construction

DNA encoding 1-169 aa (Q9NZZ3) were fused with 6His tag.

Organism

Human

Expression Host

E. coli

Activity

Not tested.

Purity

>90%, by SDS-PAGE with Coomassie Brilliant Blue staining.

Endotoxin

Please contact the lab for more information.

Stability

Samples are stable for up to twelve months from

date of receipt at -70°C

Molecular Mass

Predictes a molecular mass of 21.07 kDa. In SDS-PAGE under reducing conditions, it migrates as an approximately 24 kDa band.

Formulation

Lyophilized from sterile PBS, pH 7.4

Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Reconstitute at 0.25 µg/µl in sterile water for short-term storage. Reconstitution with 50% glycerol solution is recommended for longer term storage (see Stability and Storage for more details). If a different concentration is needed for your purposes please adjust the reconstitution volume as required (please note: the ion concentration of the final solution will vary according to the volume used). Note: Centrifuge vial before opening. When reconstituting, gently pipet and wash down the sides of the vial to ensure full recovery of the protein into solution.