

PHKA1 antibody



Catalog Number: 113818

Product name

PHKA1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PHKA1 Rabbit Polyclonal antibody. Positive WB detected in mouse skeletal muscle tissue, mouse heart tissue. Positive IF detected in HeLa cells. Observed molecular weight by Western-blot: 130 / 140 kDa

Preparation

This antibody was obtained by immunization of PHKA1 recombinant protein (Accession Number: NM_001122670). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:10-1:100

Validations



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:113818(PHKA1 Antibody) at dilution of 1:600



Immunofluorescent analysis of HeLa cells using Catalog No:113818(PHKA1 Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)