

# Anti-ICAM1 Rabbit antibody



Catalog Number: 166818

## Product name

Rabbit IgG

Anti-ICAM1 Rabbit antibody

## Applications

## Specificity

WB, IHC-P, FC, IP

Human

## Dilutions

## Antibody description

WB: 1/1000

ICAM1 Rabbit polyclonal antibody

IHC: 1/50

## Preparation

FC: 1/20-1/100

Antigen: Recombinant protein of human ICAM1

IP: 1/20

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage


Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type

## Validations



Western blot detection of ICAM1 in K562,Hela cell lysates using ICAM1 Rabbit pAb(1:1000 diluted).Predicted band size:58kDa.Observed band size:89,92kDa.