

Anti-GTF2H1 (3E8) Mouse antibody



Catalog Number: 168118

Product name

Anti-GTF2H1 (3E8) Mouse antibody

Specificity

Human

Antibody description

GTF2H1 (3E8) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant human GTF2H1 protein fragments expressed in E.coli.

Formulation

Purified mouse monoclonal antibody in buffer containing 0.1M Tris-Glycine(pH 7.4 150 mM NaCl)with 0.02% sodium azide 0.1mg/mlBSA and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications


WB, IP, ICC/IF

Dilutions


WB: 1/300

;ICC: 1/50


Validations



Immunoprecipitation analysis of HeLa cell lysates using GTF2H1 mouse mAb.



Western blot detection of GTF2H1 in HeLa, HeLa NE and Jurkat cell lysates using GTF2H1 mouse mAb (1:300 diluted). Predicted band size: 62KDa. Observed band size: 62KDa.



Immunocytochemistry staining of HeLa cells fixed with 4% Paraformaldehyde and using GTF2H1 mouse mAb (dilution 1:50).