

AXIN2 antibody



Catalog Number: 108518

Product name

AXIN2 antibody

Immunogen

[Human AXIN2 Recombinant protein \(GST tag & His tag\)](#)

Specificity

human, mouse, rat; other species not tested.

Antibody description

AXIN2 Rabbit Polyclonal antibody. Positive WB detected in human placenta tissue, mouse lung tissue. Observed molecular weight by Western-blot: 93-100kd

Preparation

This antibody was obtained by immunization of AXIN2 recombinant protein (Accession Number: BC006295). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



human placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:117079(AXIN2 antibody) at dilution of 1:200