

PC6 antibody



Catalog Number: 113618

Product name

PC6 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PC6 Rabbit Polyclonal antibody. Positive IHC detected in human brain tissue, human heart tissue, human kidney tissue, human placenta tissue, human skin tissue, human spleen tissue, human testis tissue. Positive WB detected in mouse brain tissue.

Preparation

This antibody was obtained by immunization of PC6 recombinant protein (Accession Number: NM_006200). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200


Validations



Immunohistochemical of paraffin-embedded human brain using Catalog No:113618(PC5 antibody) at dilution of 1:200 (under 10x lens)



Immunohistochemical of paraffin-embedded human brain using Catalog No:113618(PC5 antibody) at dilution of 1:200 (under 40x lens)



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:113618(PC5 Antibody) at dilution of 1:300