

Anti-IKZF3 Rabbit antibody



Catalog Number: 168028

Product name

Rabbit IgG

Anti-IKZF3 Rabbit antibody

Applications

Specificity

WB, IHC-P, ICC/IF

Human, Mouse, Rat

Dilutions

Antibody description

WB: 1/1000

IKZF3 Rabbit monoclonal antibody

IHC: 1/20

Preparation

ICC/IF: 1/50

Antigen: A synthetic peptide of human IKZF3

Validations

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA


Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type



Western blot detection of IKZF3 in HL-60,U251 cell lysates using IKZF3 Rabbit mAb(1:1000 diluted).Predicted band size:58kDa.Observed band size:58kDa.