

DDX39 antibody



Catalog Number: 109828

Product name

DDX39 antibody

Immunogen

[Human DDX39 Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

DDX39 Rabbit Polyclonal antibody. Positive IP detected in HEK-293 cells. Positive WB detected in HEK-293 cells, A431 cells, HeLa cells, HepG2 cells, Jurkat cells, mouse testis tissue, Raji cells. Positive IHC detected in human testis tissue, human kidney tissue. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 49-55 kDa

Preparation

This antibody was obtained by immunization of DDX39 recombinant protein (Accession Number: BC032128). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC, IF, IP

Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:500-1:5000


IHC: 1:20-1:200

IF: 1:10-1:100

Validations




HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:109828(DDX39 antibody) at dilution of 1:1000




Immunohistochemical of paraffin-embedded human testis using Catalog No:109828(DDX39 antibody) at dilution of 1:50 (under 10x lens)

DDX39 antibody


Catalog Number: 109828



Immunohistochemical of paraffin-embedded human testis using Catalog No:109828(DDX39 antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of HepG2 cells, using DDX39 antibody Catalog No:109828 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit



IP Result of anti-DDX39 (IP:Catalog No:109828, 4ug; Detection:Catalog No:109828 1:1000) with HEK-293 cells lysate 3000ug.