

Fascin antibody



Catalog Number: 110528

Product name

Fascin antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

Fascin Rabbit Polyclonal antibody. Positive WB detected in SH-SY5Y cells, HeLa cells, K-562 cells, mouse thymus tissue. Positive IHC detected in human tonsillitis tissue. Observed molecular weight by Western-blot: 50-55 kDa

Preparation

This antibody was obtained by immunization of Fascin recombinant protein (Accession Number: NM_003088). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations




Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:110528(Fascin Antibody) at dilution of 1:200



(under 40x lens)

SH-SY5Y cells were subjected to SDS PAGE followed by western blot with Catalog No:110528(Fascin Antibody) at dilution of 1:600



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:110528(Fascin Antibody) at dilution of 1:200 (under 10x lens)