

HUS1B antibody



Catalog Number: 111728

Product name

HUS1B antibody

Specificity

Human, Rat; other species not tested.

Antibody description

HUS1B Rabbit Polyclonal antibody. Positive WB detected in RAW 264.7 cells, human stomach tissue, PC-3 cells. Observed molecular weight by Western-blot: 31-35kd

Preparation

This antibody was obtained by immunization of HUS1B recombinant protein (Accession Number: NM_148959). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



RAW264.7 cells were subjected to SDS PAGE followed by western blot with Catalog No:111728(HUS1B antibody) at dilution of 1:500