

# Monomethylauristatin F (MMAF) Monoclonal Antibody



Catalog Number: 176728

## Product name

Monomethylauristatin F (MMAF) Monoclonal Antibody

## Specificity

Chemical

## Antibody description

Product has been tested in direct ELISA.

## Preparation

MMAF conjugated with BSA. Protein G affinity purified

## Formulation

PBS(pH7.4), lyophilized. **Reconstitution: Adding sterile water, prepare a stock solution of 2.0 mg/ml.**

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Monoclonal

## Ig Type

Mouse IgG2b,  $\kappa$


## Applications

ELISA

## Dilutions


ELISA: 1:5,000-1:20,000

## Validations



ELISA Binding Assay of Anti-MMAF Antibody to different proteins.

ELISA analysis of Anti-MMAF Monoclonal Antibody with MMAF-BSA, MMAF-KLH, BSA and KLH. Anti-MMAF Monoclonal Antibody can bind to MMAF-BSA or MMAF-KLH, but not carrier



proteins.

ELISA Binding Assay of Anti-MMAF Antibody

Immobilized MMAF ADC at 2  $\mu$ g/mL (100  $\mu$ L/well) can bind anti-MMAF mouse monoclonal antibody.