

# Anti-TLK2/PKU-ALPHA antibody



Catalog Number: 102038

## Product name

Anti-TLK2/PKU-ALPHA antibody

## Immunogen

[Human TLK2/PKU-ALPHA \(aa 397-772, His & GST Tag\) recombinant protein](#)

## Specificity

Human TLK2 / PKU-ALPHA

**Has cross-reactivity** in ELISA with Human CD64

**No cross-reactivity** in ELISA with Human CD16a; Human CD16b; Human CD32a; Human CD32b; Mouse FCGR3

## Antibody description

Rabbit polyclonal to TLK2/PKU-ALPHA

## Preparation

Produced in rabbits immunized with purified, recombinant Human TLK2 / PKU-ALPHA (rh TLK2 / PKU-ALPHA; Q86UE8-1; Leu397-Asn772). TLK2 / PKU-ALPHA specific IgG was purified by Human TLK2 / PKU-ALPHA affinity chromatography.

## Formulation

0.2 µm filtered solution in PBS with 5% trehalose

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

IHC-P

## Dilutions

IHC-P: 0.1-1 µg/ml

## Validations



TLK2 / PKU-ALPHA Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry

Immunochemical staining of human TLK2 in human testis with rabbit polyclonal antibody (0.5 µg/mL, formalin-fixed paraffin embedded sections).