

Product name

NCKAP1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

NCKAP1 Rabbit Polyclonal antibody. Positive IHC detected in human gliomas tissue, human heart tissue. Positive WB detected in mouse brain tissue, HepG2 cells, human brain tissue. Observed molecular weight by Western-blot: 115-125 kDa

Preparation

This antibody was obtained by immunization of NCKAP1 recombinant protein (Accession Number: NM_013436). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:113038(NAP1;NCKAP1 antibody) at dilution of 1:1000



Immunohistochemistry of paraffin-embedded human heart tissue slide using Catalog No:113038(NAP1;NCKAP1 Antibody) at dilution of 1:50 (under 40x lens)



Immunohistochemistry of paraffin-embedded human gliomas tissue slide using Catalog

NCKAP1 antibody



Catalog Number: 113038

No:113038(NAP1;NCKAP1 Antibody) at dilution of 1:50 (under 40x lens)