

Anti-c-Fos (7D6) Mouse antibody



Catalog Number: 168238

Product name

Anti-c-Fos (7D6) Mouse antibody

Specificity

Human, Rat, Monkey

Antibody description

c-Fos (7D6) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant human FOS protein fragments expressed in E.coli.

Formulation

Purified mouse monoclonal in PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG1


Applications

WB


Dilutions

WB: 1/500;Other 1/500

Validations




Western blot detection of c-Fos in HCT116,COS7 and C6 cell lysates using c-Fos mouse mAb (1:500 diluted).Predicted band size:62KDa.Observed



band size:62KDa.

Western blot detection of c-Fos in HeLa(TSA treated) and HeLa cell lysates using c-Fos mouse mAb (1:500 diluted).Predicted band size:62KDa.Observed band size:62KDa.



Western blot analysis of extracts from HeLa cells,untreated,50 u03bcM TSA-treated or 50 u03bcM TPA-treated for 24h,using C-FOS antibody (1:500 diluted)(upper) or GAPDH antibody (1:10000 diluted) (lower).