

RGN/SMP30 antibody



Catalog Number: 114638

Product name

RGN/SMP30 antibody

Recommended Dilution:

WB: 1:200-1:2000

Specificity

Human, Mouse, Rat; other species not tested.

IP: 1:500-1:5000

IHC: 1:20-1:200

Antibody description

RGN/SMP30 Rabbit Polyclonal antibody. Positive IP detected in mouse liver tissue. Positive WB detected in mouse liver tissue, mouse small intestine tissue, rat liver tissue. Positive IHC detected in human liver tissue, human hepat°Cirrhosis tissue, human kidney tissue. Observed molecular weight by Western-blot: 33 kDa

Preparation

This antibody was obtained by immunization of RGN/SMP30 recombinant protein (Accession Number: NM_004683). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type


Rabbit IgG

Applications


ELISA, WB, IHC, IP

Dilutions


Validations




mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:114638(RGN,SMP30 antibody) at dilution of



IP Result of anti-RGN,SMP30 (IP:Catalog No:114638, 4ug; Detection:Catalog No:114638 1:1000) with mouse liver tissue lysate 4000ug.



Immunohistochemistry of paraffin-embedded human liver slide using Catalog No:114638(RGN,SMP30 Antibody) at dilution of 1:50



Immunohistochemistry of paraffin-embedded human liver slide using Catalog No:114638(RGN,SMP30 Antibody) at dilution of 1:50