

IL33 antibody



Catalog Number: 111738

Product name

IL33 antibody

Specificity

Human; other species not tested.

Antibody description

IL33 Rabbit Polyclonal antibody. Positive IHC detected in human tonsillitis tissue. Positive WB detected in Transfected HEK-293 cells, HeLa cells, Recombinant protein.

Preparation

This antibody was obtained by immunization of IL33 recombinant protein (Accession Number: NM_001314045). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


Validations



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:111738(IL33 Antibody) at dilution of 1:700



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:111738(IL33 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:111738(IL33 Antibody) at dilution of 1:50 (under 40x lens)