

NUDT6 Polyclonal Antibody



Catalog Number: 160048

Product name

NUDT6 Polyclonal Antibody

WB 1:500 - 1:2000

IHC 1:50 - 1:200

IF 1:50 - 1:200

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to NUDT6

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-316 of human NUDT6 (NP_009014.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type


Rabbit IgG

Applications

WB, IHC, IF


Dilutions

Validations



Western blot - NUDT6 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using NUDT6 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.



Immunofluorescence - NUDT6 Polyclonal Antibody

Immunofluorescence analysis of HeLa cells using NUDT6 antibody. Blue: DAPI for nuclear staining.