

THOC5 antibody



Catalog Number: 116048

Product name

THOC5 antibody

WB: 1:200-1:2000

IP: 1:200-1:1000

Specificity

Human, Mouse, Rat; other species not tested.

IHC: 1:20-1:200

IF: 1:10-1:100

Antibody description

THOC5 Rabbit Polyclonal antibody. Positive IHC detected in human breast cancer tissue. Positive IF detected in HepG2 cells. Positive FC detected in HEK-293 cells. Positive IP detected in HEK-293 cells. Positive WB detected in HEK-293 cells, HeLa cells, NIH/3T3 cells. Observed molecular weight by Western-blot: 79 kDa

Preparation

This antibody was obtained by immunization of THOC5 recombinant protein (Accession Number: NM_001002877). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, IHC, WB, IF, IP, FC


Dilutions

Recommended Dilution:

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:116048(THOC5 Antibody) at dilution of 1:600




IP Result of anti-THOC5 (IP:Catalog No:116048, 3ug; Detection:Catalog No:116048 1:300) with HEK-293 cells lysate 2000ug.


THOC5 antibody




Catalog Number: 116048



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using Catalog No:116048(THoc5 Antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of HepG2 cells using Catalog No:116048(THoc5 Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)



1X10⁶ HEK-293 cells were stained with 0.2ug THoc5 antibody (Catalog No:116048, red) and control antibody (blue). Fixed with 4% PFA blocked with 3% BSA (30 min). Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) with dilution 1:1500.