

Anti-E2F1 (4G8) Mouse antibody



Catalog Number: 168448

Product name

Anti-E2F1 (4G8) Mouse antibody

Specificity

Human, Rat

Antibody description

E2F1 (4G8) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant human E2F-1 protein fragments expressed in E.coli.

Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG2b

Applications


WB, ICC/IF, IP

Dilutions

WB: 1/500

;ICC: 1/100


Validations



Immunofluorescent analysis of HeLa cells fixed with 4% Paraformaldehyde and using anti-E2F-1 mouse mAb (dilution 1:100). DAPI was used to

stain nucleus(blue).

Immunoprecipitation analysis of HeLa cell lysates



using E2F-1 mouse mAb.

Western blot detection of E2F-1 in C6,Raw264.7 and HeLa cell lysates using E2F-1 mouse mAb (1:500 diluted).Predicted band size:70KDa.Observed band size:70KDa.

