

Human ST6GALNAC2/SIAT7B/SIATL1/STHM (His tag) recombinant protein



Catalog Number: 518448

General Information

Gene Name Synonym

Alpha-N-Acetylgalactosaminide Alpha-2,6-Sialyltransferase 2, GalNAc Alpha-2,6-Sialyltransferase II, ST6GalNAc II, ST6GalNAcII, SThM, Sialyltransferase 7B, SIAT7-B, ST6GALNAC2, SIAT7B, SIATL1, STHM

Protein Construction

A DNA sequence encoding the Human ST6GALNAC2/SIAT7B/SIATL1/STHM Ser29-Arg374 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM PB, 150mM NaCl, pH7.4.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human ST6GALNAC2/SIAT7B/SIATL1/STHM (His tag) recombinant protein