

RXRG antibody



Catalog Number: 114948

Product name

RXRG antibody

Immunogen

[Human RXRG Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

RXRG Rabbit Polyclonal antibody. Positive WB detected in mouse heart tissue, Jurkat cells, mouse ovary tissue. Observed molecular weight by Western-blot: 51-55 kDa

Preparation

This antibody was obtained by immunization of RXRG recombinant protein (Accession Number: NM_006917). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:114948(RXRG Antibody) at dilution of 1:600