

# PEF1 antibody



Catalog Number: 113648

## Product name

PEF1 antibody

## Immunogen

[Human PEF1 Recombinant protein \(GST tag & His tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

PEF1 Rabbit Polyclonal antibody. Positive IP detected in HEK-293 cells. Positive WB detected in HEK-293 cells. Observed molecular weight by Western-blot: 30-33 kDa

## Preparation

This antibody was obtained by immunization of PEF1 recombinant protein (Accession Number: NM\_012392). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

## Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:113648(PEF1 antibody) at dilution of 1:1000



IP Result of anti-PEF1 (IP:Catalog No:113648, 4ug; Detection:Catalog No:113648 1:500) with HEK-293 cells lysate 3200ug.