

NIN Polyclonal Antibody



Catalog Number: 165748

Product name

NIN Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to NIN

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-287 of human NIN (NP_891991.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF


Dilutions

WB 1:500 - 1:1000

IHC 1:50 - 1:100


IF 1:50 - 1:100

Validations



Western blot - NIN Polyclonal Antibody


Western blot analysis of extracts of Jurkat cells, using NIN antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure



time: 90s.

Immunohistochemistry - NIN Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human stomach using NIN antibody at dilution of



1:100 (40x lens).


Immunohistochemistry - NIN Polyclonal Antibody

NIN Polyclonal Antibody



Catalog Number: 165748

Immunohistochemistry of paraffin-embedded rat kidney using NIN antibody at dilution of 1:100




(40x lens).

Immunohistochemistry - NIN Polyclonal Antibody

Immunohistochemistry of paraffin-embedded

mouse lung using NIN antibody at dilution of



1:100 (40x lens).

Immunofluorescence - NIN Polyclonal Antibody

Immunofluorescence analysis of U2OS cells using NIN antibody at dilution of 1:100. Blue: DAPI for nuclear staining.