

NDRG4 antibody



Catalog Number: 113058

Product name

NDRG4 antibody

IHC: 1:20-1:200

IF: 1:10-1:100


Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

NDRG4 Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive IHC detected in human medulloblastoma. Positive WB detected in human brain tissue. Observed molecular weight by Western-blot: 37kDa, 38 kDa and 41kDa

Validations




Preparation

This antibody was obtained by immunization of NDRG4 recombinant protein (Accession Number: NM_001242836). Purification method: Antigen affinity purified.

human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:113058(NDRG4 antibody) at dilution of 1:400

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Immunohistochemical of paraffin-embedded human medulloblastoma using Catalog No:113058(NDRG4 antibody) at dilution of 1:50 (under 10x lens)

Ig Type

Rabbit IgG


Applications

ELISA, WB, IHC, IF

Dilutions

Recommended Dilution:

WB: 1:500-1:5000



Immunofluorescent analysis of HepG2 cells, using

NDRG4 antibody



Catalog Number: 113058

NDRG4 antibody Catalog No:113058 at 1:25

dilution and Rhodamine-labeled goat anti-rabbit
IgG (red).