

Anti-NKX3-1 antibody



Catalog Number: 175058

Product name

Anti-NKX3-1 antibody

Specificity

Human

Antibody description

Mouse monoclonal antibody to NKX3-1

Preparation

This antigen of this antibody was recombinant protein

Formulation

Liquid, 1*TBS (pH7.4), 1%BSA, Preservative: 0.05% Sodium Azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Monoclonal

Ig Type

Mouse IgG2b

Applications

WB, IHC-P, FC

Dilutions

WB: 1:500-1:1,000

IHC-P: 1:50-1:200

FC: 1:100-1:200

Validations





Fig1: Western blot analysis of NKX3A on LNCaP cell lysate using anti-NKX3A antibody at 1/1,000



dilution.

Fig2: Immunohistochemical analysis of paraffin-embedded human liver tissue using anti-NKX3A antibody. Counter stained with hematoxylin.





Fig3: Immunohistochemical analysis of paraffin-embedded human prostate tissue using anti-NKX3A antibody. Counter stained with

Anti-NKX3-1 antibody




Catalog Number: 175058



hematoxylin.

Fig4: Immunohistochemical analysis of paraffin-embedded human prostate tissue using anti-NKX3A antibody. Counter stained with



hematoxylin.

Fig5: Flow cytometric analysis of PC-3 cells with NKX3A antibody at 1/100 dilution (green) compared with an unlabelled control (cells without incubation with primary antibody; purple).