

Anti-ATF6 Rabbit antibody



Catalog Number: 168058

Product name

Rabbit IgG

Anti-ATF6 Rabbit antibody

Applications

Specificity

WB, IHC-P, IP

Human

Dilutions

Antibody description

WB: 1/1000

ATF6 Rabbit monoclonal antibody

IHC: 1/100

Preparation

IP: 1/20

Antigen: Recombinant protein of human ATF6

Validations

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA


Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type



Western blot detection of ATF6 in A549,HL-60,U2OS cell lysates using ATF6 Rabbit mAb(1:1000 diluted).Predicted band size:75kDa.Observed band size:90-100kDa.