

CD16 antibody



Catalog Number: 109058

Product name

CD16 antibody

Specificity

Human; other species not tested.

Antibody description

CD16 Rabbit Polyclonal antibody. Positive IHC detected in human liver tissue, human tonsillitis tissue.

Preparation

This antibody was obtained by immunization of CD16 recombinant protein (Accession Number: NM_000569). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:50-1:500

Validations



Immunohistochemistry of paraffin-embedded human liver tissue slide using Catalog No:109058(CD16 Antibody) at dilution of 1:200 (under 10x lens).



Immunohistochemistry of paraffin-embedded human liver tissue slide using Catalog No:109058(CD16 Antibody) at dilution of 1:200 (under 40x lens).