

OAS1 antibody



Catalog Number: 113458

Product name

OAS1 antibody

IHC: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

OAS1 Rabbit Polyclonal antibody. Positive IHC detected in human liver tissue. Positive WB detected in HeLa cells, mouse brain tissue. Observed molecular weight by Western-blot: 48 kDa

Preparation

This antibody was obtained by immunization of OAS1 recombinant protein (Accession Number: NM_002534). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:113458(OAS1 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human liver using Catalog No:113458(OAS1 antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human liver using Catalog No:113458(OAS1 antibody) at dilution of 1:100 (under 40x lens)