

Product name

SNRBP antibody

Specificity

Human; other species not tested.

Antibody description

SNRBP Rabbit Polyclonal antibody. Positive IHC detected in human cervical cancer tissue, human liver cancer tissue. Positive WB detected in HepG2 cells, Jurkat cells. Observed molecular weight by Western-blot: 25-30 kDa

Preparation

This antibody was obtained by immunization of SNRBP recombinant protein (Accession Number: BC003530). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:115458(SNRBP antibody) at dilution of 1:500



Immunohistochemistry of paraffin-embedded human cervical cancer tissue slide using Catalog No:115458(SNRBP Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human cervical cancer tissue slide using Catalog No:115458(SNRBP Antibody) at dilution of 1:50

SNRPB antibody



Catalog Number: 115458

(under 40x lens)