

Anti-UCH37 Rabbit antibody



Catalog Number: 168958

Product name

Anti-UCH37 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

UCH37 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human UCH37

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, ICC/IF, IP


Dilutions

WB: 1/2000-1/10000


ICC/IF: 1/20

IP: 1/20

Validations



Immunocytochemistry of UCH37(green) in HeLa cells using UCH37 Rabbit mAb at dilution 1/50,



and DAPI(blue)

Western blot detection of UCH37 in 3T3, HeLa cell lysates using UCH37 Rabbit mAb(1:1000 diluted). Predicted band size: 38kDa. Observed band size: 38kDa.