

NNMT antibody



Catalog Number: 113258

Product name

NNMT antibody

IP: 1:200-1:1000

IF: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

NNMT Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive IP detected in HepG2 cells. Positive WB detected in mouse liver tissue, HepG2 cells, L02 cells. Observed molecular weight by Western-blot: 30 kDa

Preparation

This antibody was obtained by immunization of NNMT recombinant protein (Accession Number: NM_001372046). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IP, WB, IF

Dilutions


Recommended Dilution:

WB: 1:200-1:2000


Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:113258(NNMT Antibody) at dilution of 1:600



IP Result of anti-NNMT (IP:Catalog No:113258, 4ug; Detection:Catalog No:113258 1:300) with HepG2 cells lysate 2400ug.



Immunofluorescent analysis of (10% Formaldehyde) fixed HepG2 cells using Catalog

NNMT antibody



Catalog Number: 113258

No:113258(NNMT Antibody) at dilution of 1:50

and Alexa Fluor 488-conjugated AffiniPure Goat
Anti-Rabbit IgG(H+L)