

LIG3 Polyclonal Antibody



Catalog Number: 161658

Product name

LIG3 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to LIG3

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 730-1009 of human LIG3 (NP_039269.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:200

Validations




Western blot - LIG3 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using LIG3 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.



Immunohistochemistry - LIG3 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat heart using LIG3 antibody at dilution of 1:100 (40x



lens).


Immunohistochemistry - LIG3 Polyclonal Antibody

LIG3 Polyclonal Antibody



Catalog Number: 161658

Immunohistochemistry of paraffin-embedded human oophoroma using LIG3 antibody at dilution




of 1:100 (40x lens).

Immunohistochemistry - LIG3 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human esophageal cancer using LIG3 antibody at

dilution of 1:100 (40x lens).



Immunohistochemistry - LIG3 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse heart using LIG3 antibody at dilution of 1:100 (40x lens).