

# PYGB antibody



Catalog Number: 114358

## Product name

PYGB antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

PYGB Rabbit Polyclonal antibody. Positive WB detected in A549 cells, HeLa cells, MCF7 cells. Positive IF detected in A549 cells. Observed molecular weight by Western-blot: 97-105kd

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_002862). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:20-1:200

## Validations



A549 cells were subjected to SDS PAGE followed by western blot with Catalog No:114358(PYGB antibody) at dilution of 1:500



Immunofluorescent analysis of A549 cells using Catalog No:114358(PYGB Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)