

DCTD antibody



Catalog Number: 109758

Product name

DCTD antibody

Immunogen

[Human DCTD Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

DCTD Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive WB detected in human brain tissue, PC-3 cells. Observed molecular weight by Western-blot: 20kd

Preparation

This antibody was obtained by immunization of DCTD recombinant protein (Accession Number: XM_017007820). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:109758(DCTD antibody) at dilution of 1:500



Immunofluorescent analysis of HepG2 cells, using DCTD antibody Catalog No:109758 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).