

NDUFAF4 antibody



Catalog Number: 113068

Product name

NDUFAF4 antibody

Immunogen

[Human NDUFAF4 Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

NDUFAF4 Rabbit Polyclonal antibody. Positive WB detected in K-562 cells, MDA-MB-453s cells. Observed molecular weight by Western-blot: 21 kDa

Preparation

This antibody was obtained by immunization of NDUFAF4 recombinant protein (Accession Number: NM_014165). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:113068(NDUFAF4 Antibody) at dilution of 1:1200