

Anti-Cytokeratin 8 Rabbit antibody



Catalog Number: 167068

Product name

Anti-Cytokeratin 8 Rabbit antibody

Specificity

Human, Mouse

Antibody description

Cytokeratin 8 Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human Cytokeratin 8

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, FC, IP

Dilutions

WB: 1/5000-1/20000


IHC: 1/20

ICC/IF: 1/20-1/100


FC: 1/20

IP: 1/20

Validations



Immunocytochemistry of Cytokeratin 8(green) in HeLa cells using Cytokeratin 8 Rabbit pAb at



dilution 1/50, and DAPI(blue)

Western blot detection of Keratin 8 in MCF7, HeLa cell lysates using Keratin 8 Rabbit pAb (1:1000 diluted). Predicted band size: 54 kDa. Observed band size: 54 kDa.