

UBIAD1 Polyclonal Antibody



Catalog Number: 165468

Product name

UBIAD1 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to UBIAD1

Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human UBIAD1 (NP_037451.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB, IHC, IF

Dilutions


WB 1:500 - 1:2000
IHC 1:100 - 1:200
IF 1:50 - 1:200

Validations



Western blot - UBIAD1 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using UBIAD1 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure



time: 60s.


Immunohistochemistry - UBIAD1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human colon carcinoma using UBIAD1 antibody at dilution of 1:100 (40x lens).

UBIAD1 Polyclonal Antibody



Catalog Number: 165468




Immunohistochemistry - UBIAD1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human stomach using UBIAD1 antibody at dilution

of 1:100 (40x lens).

Immunofluorescence - UBIAD1 Polyclonal Antibody

Immunofluorescence analysis of U2OS cells using UBIAD1 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.



Immunohistochemistry - UBIAD1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human prostate using UBIAD1 antibody at dilution

of 1:100 (40x lens).

