

Anti-Tec Rabbit antibody



Catalog Number: 167168

Product name

Anti-Tec Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

Tec Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human Tec

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, ICC/IF, FC


Dilutions

WB: 1/2000


ICC/IF: 1/20-1/100

FC: 1/20

Validations



Immunocytochemistry of Tec(green) in HeLa cells using Tec Rabbit pAb at dilution 1/50, and



Western blot detection of Tec in K562 cell lysates using Tec Rabbit pAb(1:1000 diluted). Predicted band size: 74kDa. Observed band size: 58-75kDa.