

Pancreatic Lipase antibody



Catalog Number: 113568

Product name

Pancreatic Lipase antibody

ELISA, WB, IHC

Specificity

Human, Mouse, Rat; other species not tested.

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

IHC: 1:20-1:200

Antibody description

Pancreatic Lipase Rabbit Polyclonal antibody. Positive WB detected in human pancreas tissue, mouse pancreas tissue. Positive IHC detected in human pancreas cancer tissue, human pancreas tissue. Observed molecular weight by Western-blot: 52 kDa

Preparation

This antibody was obtained by immunization of Pancreatic Lipase recombinant protein (Accession Number: NM_000936). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality


Polyclonal

Ig Type


Rabbit IgG

Applications

Validations



human pancreas tissue were subjected to SDS PAGE followed by western blot with Catalog No:113568(PNLIP antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:113568(PNLIP antibody) at dilution of 1:50 (under 40x lens)