

PHF21B Polyclonal Antibody



Catalog Number: 165568

Product name

PHF21B Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to PHF21B

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 322-531 of human PHF21B (NP_612424.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF


Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200


IF 1:50 - 1:200

Validations




Western blot - PHF21B Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using PHF21B antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.



Immunohistochemistry - PHF21B Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human kidney using PHF21B antibody at dilution



of 1:100 (40x lens).

PHF21B Polyclonal Antibody



Catalog Number: 165568

Immunofluorescence - PHF21B Polyclonal Antibody

Immunofluorescence analysis of U2OS cells using PHF21B antibody at dilution of 1:100. Blue: DAPI for nuclear staining.