

HBG1 antibody



Catalog Number: 111268

Product name

HBG1 antibody

Immunogen

[Human HBG1 Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

HBG1 Rabbit Polyclonal antibody. Positive IP detected in K-562 cells. Positive WB detected in human placenta tissue. Observed molecular weight by Western-blot: 14-16 kDa

Preparation

This antibody was obtained by immunization of HBG1 recombinant protein (Accession Number: NM_000559). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:1000

Validations



human placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:111268(HBG1 Antibody) at dilution of 1:600



IP Result of anti-HBG1 (IP:Catalog No:111268, 3ug; Detection:Catalog No:111268 1:300) with K-562 cells lysate 2000ug.